

Flow Behind a Mitral Valve



Application Note PIV-003

The application focuses on the flow in the ventricle which is typically a triangular cavity (see figure below) and the PIV measurements presented here were performed in the mid section (mid depth) of this cavity. This section was chosen because of the size and shape of the ventricular cavity. The average dimensions for this cavity are $90 \times 60 \times 30$ (height \times width \times depth) resulting in a probable average two-dimensional flow (the depth playing here a minor role). The activation system controls and records a set of parameters enabling a high level of reproducibility and comparisons with clinical observations.


The flow rate through the mitral valve was measured with a Carolina Medical Model SR670 magnetic flow meter, and the atrial, ventricular and aortic pressures measured with Millar Instruments Model MPC500 piezoelectric pressure transducers. The blood flow circuit was seeded with small ($30 \mu\text{m}$ in diameter) neutrally buoyant polyamide particles with high refractive index.

The instantaneous velocity fields were determined by means of a PIV system from TSI Incorporated. This system uses a LASERPULSE™ synchronizer and PIVCAM 10-30 camera. The camera working in the cross-correlation mode was equipped with a 60mm Nikon lens. System control and data analysis was done using INSIGHT™ software.





Courtesy: Y. Knapp, E. Bertrand, F. Mouret and J. Stefanini
 Laboratoire d'Hydrodynamique Complexe, Université d'Avignon,
 Ecole Supérieure de Mécanique de Marseille



Reference

Knapp, Y., Bertrand, E., Mouret F., Stefanini, J., “Flow structure characterization behind mitral valves in a double activated heart simulator,” 5th International Symposium on Particle Image Velocimetry, Busan, Korea, September 22–24, 2003.

TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

USA	Tel: +1 800 874 2811	E-mail: fluid@tsi.com	Website: www.tsi.com
UK	Tel: +44 149 4 459200	E-mail: tsiuk@tsi.com	Website: www.tsiinc.co.uk
France	Tel: +33 491 11 87 64	E-mail: tsifrance@tsi.com	Website: www.tsiinc.fr
Germany	Tel: +49 241 523030	E-mail: tsigmbh@tsi.com	Website: www.tsiinc.de
India	Tel: +91 80 41132470	E-mail: tsi-india@tsi.com	
China	Tel: +86 10 8260 1595	E-mail: sibeijing@tsi.com	
Singapore	Tel: +65 6595 6388	E-mail: tsi-singapore@tsi.com	



TRUST. SCIENCE. INNOVATION.

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.