

Condensation Particle Counter

Model 3781

Standard replacement parts for Model 3781 Condensation Particle Counter consist of:

<u>Part Number</u>	<u>Part Description</u>
390065	Aerosol Instrument Manager [®] Software for CPC/EAD Maintenance Kit for 3781 WCPC (includes 5 wicks, 1 inline filter, 1 inlet screen,
1031545	1 transport flow orifice, and 1 pump silencer)
1031558	Inlet Cyclone for the 378X WCPCs
1500230	Wick for 3781 (box of 10)
1031544	Fill Bottle for 3781 WCPC
1181001	3781 Transport Flow Filter
1181002	3781 Damper Orifice
1181003	3781 Flow Control Orifice
1181004	3781 Transport Flow Orifice-Screenless
1181005	3781 Inlet Screen
1500229	Pump Silencer for Model 3781
1200020	Power Supply 3781 WCPC (Desktop)
1500233	Large 4L Fill and Drain Bottles for 3781/2/5/6
1031492	Large Charcoal Filter Kit for CPCs (includes 5 large charcoal filters)
1031493	Small Charcoal Filter Kit for CPCs (included 5 small charcoal filters)

For pricing information, or to place an order you can 1) send the necessary information to answers@tsi.com or 2) contact your local TSI representative (see contact information below) and speak to one of our customer specialists. Please include your contact information (telephone number and mailing address) with all requests.

Aerosol Instrument Manager[®] is a registered trademark of TSI Incorporated.

TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

USA	Tel: +1 800 874 2811	E-mail: particle@tsi.com	Website: www.tsi.com
UK	Tel: +44 149 4 459200	E-mail: tsiuk@tsi.com	Website: www.tsiinc.co.uk
France	Tel: +33 491 95 21 90	E-mail: tsifrance@tsi.com	Website: www.tsiinc.fr
Germany	Tel: +49 241 523030	E-mail: tsigmbh@tsi.com	Website: www.tsiinc.de
India	Tel: +91 80 41132470	E-mail: tsi-india@tsi.com	
China	Tel: +86 10 8260 1595	E-mail: tsibeijing@tsi.com	

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.



TRUST. SCIENCE. INNOVATION.